



# Perspectum Portal

## Customer User Guide

PDM9003

## Table of Contents

<b>PREAMBLE .....</b>	<b>4</b>
<b>1. INTRODUCTION .....</b>	<b>4</b>
<b>2. SCOPE .....</b>	<b>4</b>
<b>TERMS AND ABBREVIATIONS .....</b>	<b>4</b>
<b>DISCLAIMER .....</b>	<b>5</b>
<b>WHAT IS PERSPECTUM PORTAL? .....</b>	<b>6</b>
<b>GETTING STARTED .....</b>	<b>7</b>
<b>HOW DO I REQUEST AN ANALYSIS? .....</b>	<b>9</b>
<b>UPLOADING A NEW STUDY .....</b>	<b>9</b>
<b>REQUEST A RE-ANALYSIS .....</b>	<b>11</b>
<b>THE PROJECT TRACKER .....</b>	<b>13</b>
<b>EXPLORING STUDIES .....</b>	<b>14</b>
<b>SUBMITTED FILES .....</b>	<b>14</b>
<b>EXPLORING ANALYSES .....</b>	<b>15</b>
<b>STUDY COMMUNICATIONS .....</b>	<b>15</b>
<b>WHAT QUALITY CHECK CAUTIONS MEAN .....</b>	<b>16</b>
<b>HOW CAN I MANAGE MY ACCOUNT? .....</b>	<b>17</b>
<b>CHANGING YOUR EMAIL ADDRESS .....</b>	<b>17</b>
<b>CHANGING YOUR PASSWORD .....</b>	<b>17</b>
<b>MANAGING YOUR PROJECTS .....</b>	<b>17</b>
<b>ACTIVITY LOG .....</b>	<b>18</b>
<b>SUBSCRIBING AND UNSUBSCRIBING TO PROJECT NOTIFICATIONS .....</b>	<b>18</b>
<b>SETTING UP NEW PROJECTS .....</b>	<b>18</b>
<b>SECURITY – KEEPING YOUR DATA SAFE .....</b>	<b>19</b>
<b>ADDING NEW USERS .....</b>	<b>19</b>

ACCOUNT ACTIVATION..... 19

CHANGING EXISTING USER PERMISSIONS..... 20

BROWSER REQUIREMENTS..... 20

I AM HAVING PROBLEMS WITH PERSPECTUM PORTAL – WHAT SHOULD I DO? ..... 21

FORGOTTEN PASSWORD..... 21

I CAN’T SEE THE DATA I AM EXPECTING TO SEE ..... 21

SOMETHING ELSE? ..... 21

## 1. Introduction

This document is intended to be used as a customer guide to using Perspectum Portal, Perspectum's online Portal for uploading data to its Quantitative Analysis Service (QAS).

## 2. Scope

This document provides general instructions for the use of the Perspectum Portal, including legal disclaimers, an overview of the key features, a guide to the workflow, information required for setup and FAQs.

## Terms and Abbreviations

Term	Meaning
<b>Blinded report</b>	A 'blinded report' is a report that users can't access without generic user download permission. This report status is used within clinical trials to limit access to trial data.
<b>Dashboard</b>	The 'homepage' of Perspectum Portal. Users can see all Studies from Projects that they have access to here
<b>MRI</b>	Magnetic Resonance Imaging
<b>Organisation</b>	A business (e.g. University, Pharmaceutical company, Radiology Centre) who own or contribute to Projects
<b>Patient Identifiable Information</b>	Information which is private to an individual and which could enable the identification of that individual. For example, name and address or information about health
<b>Perspectum Portal</b>	Formerly Edison. Perspectum's cloud-based Portal for data transfer and delivery.
<b>Perspectum Portal User</b>	Anyone who has access to Perspectum Portal. A Perspectum Portal User must be affiliated to an Organisation
<b>Project</b>	A collective term for a group of related MRI datasets, e.g. a clinical service or a pharmaceutical trial, identified as such for ease of billing and access. All Studies must belong to a single Project
<b>QAS</b>	Quantitative Analysis Service
<b>Study</b>	A dataset from a single patient and a single scanning session. Studies are uploaded to Perspectum Portal for analysis

Please read this Customer User Guide (the 'Guide') carefully before using Perspectum Portal. By continuing to use the Perspectum Portal, you accept the terms of this disclaimer.

Whilst Perspectum will use reasonable efforts to ensure the accuracy of all the information in the Guide to the extent permitted by law, Perspectum assumes no liability to any party for any damage or loss caused by errors or omissions or statements of any kind in the Guide, its updates or supplements, whether such errors are omissions or statements resulting from negligence, accidents or any other cause. Perspectum assumes no responsibility for the content of external website links or documents referred to within the Guide and the use of these should not be interpreted as an endorsement by Perspectum of those links and/or documents.

To the extent permitted by law, Perspectum assumes no liability arising from the application or use of the Guide; nor any liability for incidental or consequential damages arising from the use of the Guide. Perspectum disclaims all warranties regarding the information contained herein whether expressed or implied.

Perspectum is committed to keeping the information in the Guide up to date and correct. In the event of a material change to the Guide, Perspectum will issue you with an updated version as soon as reasonably practicable. Any reliance you place on an outdated printed version of the Guide is strictly at your own risk.

Whilst Perspectum will use its reasonable efforts to keep the website up and running smoothly we do not take responsibility for, and will not be liable for, the website being temporarily unavailable due to technical issues beyond our control.

To the extent permitted by law, Perspectum shall not be liable for any loss or damage including without limitation, indirect or consequential loss or damage, or any loss or damage whatsoever arising from loss of data or profits arising out of, or in connection with, the use of the Guide.

If you have any questions about the Guide, please contact Perspectum's support team via email at [support@perspectum.com](mailto:support@perspectum.com).

Perspectum Portal is the interface through which Perspectum's Quantitative Analysis Service (QAS), which includes *LiverMultiScan*, is delivered to our customers. Its key features are shown in Figure 1.




Figure 1 Key features of Perspectum Portal

Getting started

When you first log on to Perspectum Portal, you are taken to the **Dashboard** (Figure 2). This shows all the patient studies that have been uploaded that you have access to. Each row refers to an individual patient dataset (referred to as a ‘Study’) defined by its **Study ID**. Every Study is associated with a **Project**. ‘**Projects**’ can be clinical services, research studies or clinical trials.

Patient Studies

Project Tracker

Upload New Study

All Fields Search for... Data from all time...

Legend

Study ID ▲	Project ▲	Organisation ▲	Uploaded By ▲	Date Uploaded ▲	Analyses	Reports
Test study 3	SQA Project 1	SQA Hospital A	Customer SQA E2E	Thu 17 Feb 2022, 02:10 pm	Simple ✓ Advanced ✓	Download reports
Test study 2	SQA Project 1	SQA Hospital A	Customer SQA E2E	Thu 17 Feb 2022, 02:07 pm	Simple	-
Test study	SQA Project 1	SQA Hospital A	Customer SQA E2E	Thu 17 Feb 2022, 02:06 pm	Simple i Advanced ✓	Download report

Show 10 50 100 results

Figure 2 Dashboard page of Perspectum Portal listing all patient studies that have been submitted for analysis.

The dashboard shows whether the report is available for each Study via the ‘**Reports**’ column. Where it is available you will see the option to ‘**Download report**’.

The Dashboard can be filtered by Study ID, Project, Analysis Type, Organisation, Date Uploaded and who uploaded the study using the ‘**Search for.....**’ search bar. You can also search for datasets uploaded over specific time periods, by clicking on the adjacent ‘**Showing data from all time...**’ window (Figure 3). Both text and data searches can be exited by clicking on the cross icon at the end of ‘**Search results for:**’.

All Fields Search for... Data from all time...

Legend

This week  
Last week  
This month  
Last month  
Year to date  
Custom Range

2022-02-17

2022-02-17

< Feb 2022 Mar 2022 >

Su	Mo	Tu	We	Th	Fr	Sa	Su	Mo	Tu	We	Th	Fr	Sa
30	31	1	2	3	4	5	27	28	1	2	3	4	5
6	7	8	9	10	11	12	6	7	8	9	10	11	12
13	14	15	16	17	18	19	13	14	15	16	17	18	19
20	21	22	23	24	25	26	20	21	22	23	24	25	26
27	28	1	2	3	4	5	27	28	29	30	31	1	2
6	7	8	9	10	11	12	3	4	5	6	7	8	9

Analyses	Reports
Simple ✓ Advanced ✓	Download reports
Simple	-
Simple i Advanced ✓	Download report

Show 10 50 100 results

Figure 3 Search by date or time period

The number of studies visible at a time on the Dashboard can be altered. On the bottom right of the list of studies on the Dashboard, use the number selectors to change the number of studies that are shown, as shown in Figure 4. Currently, either 10, 50 or 100 studies can be selected for display at a time on the Dashboard.

Final •  
Restricted

PDM9003-v25.1

Page 7 of 22  
Uncontrolled if printed

Test study	SQA Project 1	SQA Hospital A	Customer SQA E2E	Thu 17 Feb 2022, 02:06 pm	<span>Simple</span>   <span>Advanced ✓</span>	<input type="button" value="Download report"/>
Show <a href="#">10</a> <a href="#">50</a> <a href="#">100</a> results						

Figure 4 Selecting the number of Studies to display on the Dashboard



## How do I request an analysis?

### Uploading a new study

1. On the **Dashboard** (Figure 2), click on the 'Upload New Study' link

### Select Project

2. At the Upload Patient Study page (Figure 5) choose the **Project** that the Study is associated with (Figure 6). This may be a specific clinical trial, research project or clinical service project

### Label your study

3. For clinical trials, where anonymised subjects are scanned at multiple timepoints, the Project will have been set up to expect both a **Patient ID** and a **Timepoint** entry for each scan. For such projects, both these fields will appear on the 'Upload' page and are mandatory. Select an existing **Patient ID** or create a new **Patient ID** following any specified naming conventions, ensuring it does not contain any Patient Identifiable Information. Select a **Timepoint** from the drop-down menu. The **Study ID** will be automatically generated based on the **Patient ID** and **Timepoint** selected.
4. Other kinds of Project (e.g. clinical scans or research projects without multiple timepoints) don't require a Patient ID and a Timepoint, rather just a single **Study ID**. This Study ID must not contain any Patient Identifiable Information; however, it may contain information that identifies the scan, such as an accession number. It could also be a deidentified patient ID.

### Select the Uploading and Referring Organisations

5. Select the **Uploading Organisation** (usually the Imaging Centre where the data was acquired) and, if required, the **Referring Organisation** (usually the clinical site where the referral for the scan has come from) that the **Study** is associated with. You will only be asked to select a Referring Organisation, if they have been set up for your selected project.
6. Select the applicable **Analysis Type** (note that more than one can be selected) (Figure 7). Note that if you forget to add an Analysis Type, it can be added later in the study view page by users with the required permissions.
7. Provide any additional details or requests for analysis in the comments window.
8. Drag and drop the patient study into the file upload window at the bottom of the page (Figure 5). We highly recommend zipping the files/folders of the patient study in to one zip file as it will upload much faster than adding the individual files.
9. Press the **Submit** button.
10. The data will now appear as a New Study in the **Dashboard**.

[Dashboard](#) / [Upload Patient Study](#)

**Project:**

CoverscanMD Test

Select the Project that the study relates to

**Uploading Organisation:**

Perspectum RND

Select the Organisation that the study was generated at

**Referring Organisation:**

N/A

Select the Organisation that referred this study

**Study ID:**

Study ID

Enter the study identifier here. **Please ensure that NO patient identifiable information is included in the Study ID**

**Comments:**

Comments

Add any comments to be recorded with the study. **Please ensure that NO patient identifiable information is included in the comment.**

**Analysis Types:**

CoverScan MD

Select one or more (hold ctrl/cmd and click) types of Analysis that this Study should undergo

To add a patient (study), drop files here, or click to select files manually

- All files should be from the same study
- Individual files or a ZIP are fine

Submit

Figure 5 Upload Patient Study page. File upload window is at the bottom of the page. Note that some of these fields may not be available for your Project.

**Project:**

✓ Select project

Test Project

Project Too

Quality Control

Study ID

Enter your ID for the study

Figure 6 Select project drop down menu

Analysis Types:

Advanced

CoverScan MD

Select one or more (hold ctrl/cmd and click) types of Analysis that this Study should undergo

Figure 7 Select analysis type

Request a Re-analysis

If you have been given permission to request a study re-analysis, and the project allows it, you can do so via the **Summary** tab (Figure 8).

Summary

Submitted Files

Comments

History

System Administration

Summary

Project

Ines Proj

Project Group

None

Uploading Organisation

Ines Org

Project Owner

Ines Org

Uploaded by

Ines Almeida (Ines.Almeida@perspectum.com)

Data Rights Position

Must consult Legal

Add analyses

Request study re-analysis

Figure 8 Reanalysis requests are made from the Summary tab

After clicking “Request study re-analysis”, a window will appear showing the re-analysis options (Figure 9). All available fields must be completed to continue with your re-analysis request. Once you have requested a Re-analysis, a new entry will appear on your dashboard which will be processed as another study.

Re-analyse Study

Please select a project from the options below:

CoverscanMD Test

Please change Study ID below, if required:

Study ID:

Study\_012

Please ensure that NO patient identifiable information is included in the Study ID.

Please select the type of re-analysis from the options below:

☐ Standard

☒ Advanced

☐ Enhanced cT1

☐ LMS

Select one or more types of Analysis

A comment is required for re-analysis:

Enter comment...

Please ensure that NO patient identifiable information is included in the comment.


Enter password...


Cancel Re-analyse Study

Figure 9 Request a re-analysis

## The Project Tracker

If a Project Tracker has been configured for your project, and you have permission to view it, you can explore it by pressing the “Project Tracker” button on the **Dashboard** (Figure 2). The Project Tracker allows you to take a bird’s eye view of your whole project. From here you can follow the progress of all patients within your project (in progress, delivered or returned without analysis), and navigate to individual analyses of interest, across multiple timepoints. You can also filter by Analysis Type, if desired. You can also export your project tracker by clicking the Export button shown in Figure 10. Note that the Analysis Type filter applies to the Export as well, so the .csv file will only contain the selected Analysis Type.





Perspectum Portal 


**Project Tracker**


Logged in as ines@perspectum.com  
(Sign out)  
Last Login on Tue 23 Mar 2021, 11:35 am  
[Manage Account](#)






Dashboard / Project Tracker

Patient ID Search for... 

Viewing Project: Timepoint Proj 

Viewing Analysis Type: - 

Export 

Patient ID ▲	Timepoints 		
	A	B	C
Patient_1	 	 	

Study 33: LMS Basic

Figure 10 The Project Tracker

## Exploring Studies

Clicking on the related **Study ID** for a study on the **Dashboard** opens the **Study** page (Figure 11). Here you can find details of the Study including the **Project**, **Project Owner**, who it was **Uploaded by**, the **Uploading Organisation**, the **Referring Organisation**, **Patient ID** (if any), **Timepoint** (if any), as well as the **History** of updates to the Study.

Each analysis is displayed on this page – there is a button to download the report (if available), as well as a button to view the analysis in more detail (see below).

Any comments attached to the Study can be accessed via the **Comments** tab. If there have been any issues with the quality of data acquisition, these will be available here. You can also filter comments by Analysis Type.




Figure 11 Overview of an individual study

## Submitted files

The data submission can be viewed in the **Submitted Files** tab (Figure 13). Perspectum Portal pulls information directly from the uploaded **Study**, including the scanner model and software version that the data was acquired with, together with the date that the **Study** was acquired and the **Patient ID** (if entered on the scanner).

Summary

Submitted Files

Comments

History

System Administration

Study Files

Download all

DICOM Study 1

Scanner Software Version	syngo MR E11
Instance Creator UID	Tag does not exist
Scanner Manufacturer	SIEMENS
Scanner Model	Aera
Patient Name	TD053_GEMINI
Patient ID	TD053_GEMINI
Magnetic Field Strength	1.5
Technician	AB
Device Serial Number	142048
Study Date	20201012

Figure 12 DICOM Study overview

Exploring Analyses

Within the Analysis Details section, you can find more details relating to each individual analysis, including Data Files, and Reports.

Dashboard / Study Test study 2 / Analysis Details

Summary

Analysis Links

Submitted Files

Comments

Data Files

Reports

Summary

Project

SQA Project 1

Project Owner

SQA Hospital A

Uploaded by

Customer SQA E2E

Uploading Organisation

SQA Hospital A

Analysis: Simple

Submitted

Under Analysis

Under Review

Report Available

Figure 13 Overview of an individual analysis

Study Communications

The **Comments** tab is where communications about a **Study** are displayed. You can also leave comments about a specific **Analysis**, if applicable.

Final •  
Restricted

PDM9003-v25.1

Page 15 of 22  
Uncontrolled if printed

We may also send an email alerting you that a response is required. An example email is shown in Figure 14.

Hello Ben,

A comment has been added to the study **Test study**, which was uploaded by you on Perspectum Portal.

Please see [https://portal.perspectum.com/view/78476?tab=comments&comment\\_id=1082004](https://portal.perspectum.com/view/78476?tab=comments&comment_id=1082004) to view the new comment.

Figure 14 Example email requesting a response from a Customer.

## What Quality Check Cautions mean

If configured for your Project, Perspectum Portal may detect certain issues with a **Study**, as shown in Table 1.

Table 1. Quality Check Cautions

Quality Issue	Probable cause	Recommended action
No DICOM Data	No DICOM data was present in the uploaded study data	Re-upload the full dataset
Prior Upload	The unique study identifier of your study has been previously uploaded to Perspectum Portal	Check that the uploaded data is correct and reupload if required
Missing Series	Incorrect or incomplete acquisition or upload	Check that series data was acquired and reupload
Duplicate Data	Study files have been duplicated within a dataset upload	Check that the uploaded data is correct and reupload if required
Multiple Studies	More than one unique MRI study uploaded in a single upload	Reupload individual studies separately

**Quality Check Cautions** are shown in the **Submitted Files** tab (e.g. Figure 15). Clicking on one of the Quality Check Cautions (e.g. Missing Series) provides further information on the problem.

Summary

**Submitted Files**

Comments

History

System Administration

### Study Files

**Quality Check Cautions**

Potential problems with the Study have been detected

QC details

Missing Series Prior Upload

One or more required tags missing

Figure 15 Caution where Missing Series, as well as a Prior Upload, has been detected



## How can I manage my account?

Your account management page can be reached by clicking on **Manage Account** at the top right of the **Dashboard**. Here you can change your email address and password as well as checking which **Projects** you have access to and recent account activity (Figure 16)




Figure 16 Account management page. Top left: Customer details pane.

### Changing your email address

To change your email address, click **Change Email Address** in the **Customer Details** pane. You will need to enter your new email address, together with your password.

### Changing your password

To change your password, click **Change Password** in the **Customer Details** pane of the **Manage Account** page. You will need to enter your existing password as well as your new password. To keep your data as safe as possible, passwords must meet these criteria:

- Contain at least 8 characters
- Contain at least one lowercase character, one uppercase character and a number
- Not be entirely numeric

Passwords should be known only to the person that they have been assigned to. They should not be written down, emailed or passed among colleagues. If you suspect that your password has been compromised, change it as soon as possible, and contact [support@perspectum.com](mailto:support@perspectum.com)

N.B. Passwords must be updated periodically to ensure ongoing security of data.

### Managing your projects

To manage your projects, click **Manage Projects** (Figure 16). This opens the list of **Projects** that you have access to, together with your individual permissions and notifications for that **Project**. Depending on your role, you may have permission for all, some or none of the following actions:

- Download Blinded Report – User can still download reports that have been set to blinded
- Download Dataset – User can download the data uploaded to the study
- Download Data Files – User can download any additional data files uploaded to the study
- Download Report – User can download reports uploaded to the study (that have not been set to blinded)
- Re-analyse Study – User can request re-analysis of a study
- Upload Dataset – User can upload data to the Portal during study upload
- View Project Tracker – User can view the Project Tracker for that project, if available

You can also check whether you are receiving notification emails:

- 'Daily': the status of one or more of your Study Datasets has changed in the last twenty-four hours (daily update).
- 'Rejection': all analyses from a study have been rejected (real-time as soon as it is rejected).
- 'Quality issues': a Study Dataset has raised one or more quality check cautions.
- 'Report Delivered': a report has been delivered on a study.

## Activity log

To view the recent activity on your account, click **View Activity** (Figure 16). In the Activity screen you can then see:

- Previous logins (successful and unsuccessful)
- Changes of password or email
- When you uploaded **Studies**
- When you viewed or downloaded **Reports** and data.

## Subscribing and unsubscribing to Project Notifications

You can change your notifications preferences for all **Projects** to which you have access at any time. Simply check or uncheck the notification boxes to the right of the **Project**.

Dashboard / Manage Account / Manage Projects			
Projects			
Project	Owning Organisation	Permissions	Notifications
OBH_eCRF_Prj	OBH_eCRF_Org	Download Report Upload Dataset	<input type="checkbox"/> Daily <input type="checkbox"/> Rejection <input type="checkbox"/> Quality Issues <input type="checkbox"/> Report Delivered

Figure 17 Project section on Account Management page. You can subscribe or unsubscribe to available projects at any time.

## Setting up new Projects

To add a new **Project** to your **Organisation**, you will need to contact a member of our Service team, who will guide you through the setup process, to ensure we have the appropriate details. For each **Project** that your **Organisation** owns, the following information is required:

1. A list of **Staff users** who have permission to view the **Project**

2. For these staff, the following **Permissions/Preferences** must also be defined
  - a. Upload Study
  - b. Download Report
  - c. Download Blinded Report
  - d. Request Re-analysis
  - e. Rejections notification email
  - f. Daily update email
  - g. Add analysis/request re-analysis
3. A list of **Organisations** (excluding yours, the owner **Organisation**) who also have access to the **Project** (partner **Organisations**).

For each **Project** that your **Organisation** contributes to but does not own, the information in points 1 and 2 of the above list will still be required. This may need to be confirmed by the **Organisation** who owns the **Project** (e.g. a pharmaceutical company sponsoring a multicentre trial).

Note: Where there are multiple sites contributing to a **Project** owned by another **Organisation**, members of the individual sites will only be able to see **Project** data uploaded at their site.

## Security – keeping your data safe

Data security is of paramount importance. All Perspectum Portal users have their own individual usernames and passwords. Permissions for overview, data upload and report download are set on a Project-by-Project basis and can be updated on request.

Perspectum Portal contains sensitive information. As such, the utmost care should be taken to avoid unauthorised access. To do this, please follow this guidance:

- Perspectum Portal should only be accessed by authorised users.
- You should avoid keeping a record (e.g. on paper, software file or hand-held device) of your password. If you believe that your password has been compromised, change it immediately and contact the Support desk.
- Passwords should be secure (see **Changing your password** for guidance).

## Adding new users

### Account activation

Perspectum Portal is accessed through <https://portal.perspectum.com>. If you haven't accessed Perspectum Portal before, you will first need to activate your account. You should have received, or will soon receive, an email similar to that shown in Figure 18. It is important that you follow the instructions to complete activation within the time given.




Figure 18 Example activation email

If you need further assistance, please contact your account manager via email [support@perspectum.com](mailto:support@perspectum.com) or telephone +44(0) 1865 655342.

Any changes to **Project** access require written confirmation from the **Organisation's** primary contact.

## Changing existing user permissions

Access to a **Project** can only be granted by the **Organisation** that owns that **Project**, or an authorized delegate. When **Projects** are initially set up, we ask the **Organisation** to provide details of **Staff Users** at contributing **Organisations**. New **Staff Users** can be added at any time by the **Project**-owning **Organisation** contacting Perspectum.

The primary contact at the site should email [support@perspectum.com](mailto:support@perspectum.com) with the user's name (first and last), email address, what **Projects** they need access to and what permissions they will require. The same procedure is required for changing access rights.

## Browser requirements

Perspectum Portal currently supports the internet browsers as defined below:

- Google Chrome on Windows and Mac Operating System (OS).
- Firefox on Windows and Mac OS.
- Microsoft Edge on Windows
- Safari on Mac OS.

Internet Explorer and Opera are not supported by Perspectum Portal. Note that browser support is subject to change, particularly when historic browsers reach end-of-life. **If you are unsure of your browser, please contact your system administrator.**

## I am having problems with Perspectum Portal – what should I do?

### Forgotten password

At the login page, click on the 'Forgotten Password?' link. This will send an email to your registered address with further instructions.

### I can't see the data I am expecting to see

To be able to see a **Study**, you must have both membership of that **Project** AND membership of an **Organisation** that is connected to that **Project**. There are four categories of **Organisational** access:

**1. Project-Owning Organisation**

If you are a member of an **Organisation** who owns a **Project** AND you are a member of that **Project**, you can see ALL data within that **Project**, even if it comes from **Organisations** outside of your own.

**2. Project-Uploading Organisation**

If you are a member of an **Organisation** who uploads **Studies** to (but does not own) a **Project** AND you are a member of that **Project**, you can see data within the **Project** that has been uploaded by YOUR **Organisation** only. You cannot see data that comes from an **Organisation** that you are not a member of.

**3. Project-Partnering Organisation (e.g. a Contract Research Organisation)**

If you are a member of an **Organisation** who is Partnering with an **Organisation Project** AND you are a member of that **Project**, you can see ALL data within that **Project**, even if it comes from **Organisations** outside of your own.


**4. Project-Referring Organisation**

If you are a member of an **Organisation** who refers studies for upload to (but does not upload or own) a **Project** AND you are a member of that **Project**, you can see data within the **Project** that has been referred by YOUR **Organisation** only. You cannot see data that comes from an **Organisation** that you are not a member of.

### Something else?

If you encounter problems whilst using Perspectum Portal, you can contact us through the **Report a problem** link at the bottom of the **Dashboard** page. This brings up a form where you can summarise the problem (Figure 19) To enable us to provide a solution as quickly as possible, please provide as much information as possible. We would recommend including:

- What you were attempting to do
- What you were expecting to happen
- What actually happened
- Any error message you saw
- What web browser you are using to access Perspectum Portal



The screenshot shows a web form titled "Report a problem". At the top, there is a blue sidebar with a white circle containing the letter "i". To the right of the sidebar, the text reads: "Please provide your feedback below, providing as much detail as you can, including:" followed by a bulleted list: "• What you were attempting to do", "• What you were expecting to happen", "• What actually happened", and "• Any error message you saw". Below this, there are two text input fields: "Problem summary\*" and "Detailed description\*", both with red asterisks indicating they are required. The "Detailed description" field is larger and has a small icon in the bottom right corner. Below the "Detailed description" field, there is a small blue icon with a white "i" and a question mark. Below this, there are two text input fields: "Your name" and "Your email". Below the "Your email" field, there is a checkbox with the text: "Include data about your current environment, like the browser and page URL. This helps us understand your feedback better." Below the checkbox, there is a link: "What is included in the data about my current environment?". At the bottom right of the form, there are two buttons: "Submit" and "Close".

Figure 19 The Report a Problem form

Alternatively, please contact [support@perspectum.com](mailto:support@perspectum.com) or by telephone on +44(0) 1865 655342.